

Product Name :
Pulveriser

Product Code :
NLAB-AGRETLAB160022



Description :

Pulveriser

Technical Specification :

Designed for grinding materials to produce fine mesh samples Ideal for use in Cement and Chemical Industries

Self contained grinder with a rotating disc having planetary movement in vertical plane

Reduces about 450g quartz type material to 100 mesh in one minute

Dia. of grinding wheel: 175mm

Material of grinding wheel: High Carbon Steel having 53-60 HRC hardness

Input of material: 6mm

Suitable for operation with 440V, Three Phase, A.C. supply

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com