


**Product Name :**  
Punch Press

**Product Code :**  
NLAB-MACHINELABWK26005



**Description :**

Punch Press

**Technical Specification :**

Punch Press .

Hand operated bench mounting standard single station punch press, a multi-purpose precision machine which perforates holes of various sizes, and shapes upto 4" dia in 16 gauge steel-rapidly and efficiently. Used also for an unlimited variety of blanking, drawing, embossing and forming operations. Adjustable side and back gauges incorporated in the machine allow precision gauging for exacting duplicating operations. Perfect alignment of the punch head is assured at all times by a dove-tail shaped ram which is hardened and precision ground. The roller cam provides tremendous pressure at the point of impact.

**Standard Equipment :**

Standard equipments includes punch holder, for punches with 1/2" dia, shank die holder, short and long handle, for operator convenience. Side and back gauges.

**Extra Equipments :**

Stand, Different types of dies & punches.

**Depth of throat**  
**Rated Capacity**  
**Height of throat**  
**Stroke**  
**Hole in Ram (Diameter)**  
**Bed Dimensions**  
**Slug Waste Hole**  
**Floor Space (on stand)**  
**Weight approx.**

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003