


Product Name :
Punch Press

Product Code :
NLAB-MACHINELABWK26005



Description :

Punch Press

Technical Specification :

Punch Press .

Hand operated bench mounting standard single station punch press, a multi-purpose precision machine which perforates holes of various sizes, and shapes upto 4" dia in 16 gauge steel-rapidly and efficiently. Used also for an unlimited variety of blanking, drawing, embissing and forming operations. Adjustable side and back gauges incorporated in the machine allow precision gauging for exacting duplicating operations. Perfect alignment of the punch head is assured at all times by a dove-tail shaped ram which is hardened and precision ground. The roller cam provides tremendous pressure at the point of impact.

Standard Equipment :

Standard equipments includes punch holder, for punches with 1/2" dia, shank die holder, short and long handle, for operator convenience. Side and back gauges.

Extra Equipments :

Stand, Different types of dies & punches.

Depth of throat
Rated Capacity
Height of throat
Stroke
Hole in Ram (Diameter)
Bed Dimensions
Slug Waste Hole
Floor Space (on stand)
Weight approx.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003