

Product Name :
Puritan Bennett Ventilator

Product Code :
NLABMEDICAL-M3390004

Description :

Puritan Bennett Ventilator

Technical Specification :

Puritan Bennett Ventilator

State-of-the-art gas delivery system.

Demonstrated track record of reliability.

Full-featured microprocessor-controlled ventilator with 7202 display screen.

Gas delivery features - Flow-by®, Pressure Support, Pressure Control Ventilation.

Patient monitoring features - Graphics 2.0, Respiratory Mechanics, Auto-PEEP.

Communications features - RS-232, Digital Communications Interface.

Optional Metabolic Monitor.

Integrated safety systems.

Optional humidifier is available with this unit.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003