

Product Name : Puritan Bennett Ventilator	Product Code : NLABMEDICAL-M3390004
Description :	
Puritan Bennett Ventilator	
Technical Specification :	
Puritan Bennett Ventilator	
State-of-the-art gas delivery system.	
Demonstrated track record of reliability.	
Full-featured microprocessor-controlled ventilator with 7202 display screen.	
Gas delivery features - Flow-by®, Pressure Support, Pressure Control Ventilation.	
Patient monitoring features - Graphics 2.0, Respiratory Mechanics, Auto-PEEP.	
Communications features - RS-232, Digital Communications Interface.	
Optional Metabolic Monitor.	
Integrated safety systems.	
Optional humidifier is available with this unit.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003