

**Product Name :**  
Pycnometer

**Product Code :**  
NLAB-TESTINGLABM109031



**Description :**

Pycnometer

**Technical Specification :**

Specific Gravity is the ratio of the weight in air of a given volume of a material at a standard temperature to the weight in air of an equal volume of distilled water at the same stated temperature. The Pycnometer Bottle can be used to test a wide range of materials from clay to sand and gravel smaller than 10mm. It consists of 1kg Glass Jar with Brass Cone, Locking Ring and Rubber Seal.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003