

Product Name :
Quantab Chloride Titrator Type

Product Code :
NLAB-AGRETLAB160026



Description :

Quantab Chloride Titrator Type

Technical Specification :

This method provides a measurement of water soluble chloride salts present in a soil or aggregate sample. It is based on the Volhard method.

Quantab chloride titrator, type 1175 (711195), range 0.005% to 0.1% NaCl. Pack of 40 strips.

Quantab chloride titrator, type 1176 (711196), range 0.05% to 1% NaCl. Pack of 40 strips.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003