

Product Name :
Radial Heat Conduction

Product Code :
DM184



Description :

Radial Heat Conduction

Technical Specification :

Radial Heat Conduction

Disk of brass with thickness of 3 mm and diameter of 110mm, provided with a central copper cylinder of 14 mm; the whole sample is coated with plastic material to minimize heat losses and to protect operators from scalds

Cartridge resistance of 100 W @ 24 Vdc with safety thermostat

Copper pipe around disk circumference for water cooling, provided with reducer 1/2" fi lter and control valve

Six K-type thermocouples arranged radially on disk

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003