

<b>Product Name :</b> Radiant Warmer with Under Surface Phototherapy	<b>Product Code :</b> NLABMEDICAL-M3440003
---	---

--

**Description :**

Radiant Warmer with Under Surface Phototherapy

**Technical Specification :**

Easy to use

Focused heating

Servo control mode

Air control mode

Manual mode

Under Surface Unit:

Detachable source unit can be used along with radiant warmer & infant trolley

Two, 2 feet TLD daylight tubes

Two, 2 feet imported low pressure mercury vapour fluorescent tubes emitting blue radiation with a peak at 450 nm specially used for hyperbilirubinaemia treatment (TL 20W/52)

M.S. fabrication with powder coated paint

Clear acrylic filter

Source cooling fan

**Naugralabequipments**  
**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

