

Product Name :
Ray Optics Kit

Product Code :
NLAB-EDUSCHOOLEQ1740115



Description :

Ray Optics Kit

Technical Specification :

A Ray box kit also known as a Light box kit is used to understand various optical properties .These optical properties are made practically visible and more erudite using a ray box/light box kit. A ray box/light box kit is a simple enclosed box (one side open) consisting of an incandescent bulb, various detachable slides, lenses, mirrors and color filters. These slides and lenses when used individually in a ray box/light box produce various optical properties which are determined by analyzing the distinct paths of light rays coming out of the ray box/light box.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003