


Product Name :
Reactivity Test Apparatus

Product Code :
NLAB-CEMETLAB130029



Description :

Reactivity Test Apparatus

Technical Specification :

Ground quicklime shall be tested for their reactivity with this apparatus which consists essentially of a Dewar vessel, 1000 ml cap., stirring motor, calibrated thermometer, stand and accessories.

Two versions are available:

supplied with digital thermometer

with advanced datalog thermometer, temperature probe, serial cable for PC connection and dedicated software for download data

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003