

**Product Name :**  
Rectangular Horizontal Autoclave

**Product Code :**  
NLAB-LABORATORYEM29005



**Description :**

Rectangular Horizontal Autoclave

**Technical Specification :**

Rectangular Horizontal Autoclave Inner chamber Lid & boiler is made of stainless steel. Sterilization cycle is controlled by three way valve or multi port valve made of stainless steel fitted with automatic pressure control switch, safety valve, vacuum valve, water level indicator also fitted in sterilizer. The body is mounted on the tubular stand. WORKING CHAMBER (W x H x D)

600x600x1200mm

600x900x1500mm Optional

Accessories:-

Water Cut Off Device.

Vacuum Breaker Valve

Musical Bell

Timer

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003