

**Product Name :**  
Reflux Extractor 4000 Gms

**Product Code :**  
NLAB-BTMENLAB150019



**Description :**

Reflux Extractor 4000 Gms

**Technical Specification :**

The simple apparatus working on the same operation principle of consisting of cylindrical glass jar supporting two metal cones of stainless steel cloth and a metal condenser on top of the jar. Supplied complete with 100 filter papers & wire gauge, Hot Plate.

Note: Spare Cylindrical glass jar can be supplied at an extra cost.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003