


**Product Name :**  
Refraction Kit with Cylindrical Light Box

**Product Code :**  
NLAB-SCHOOLABS20005



**Description :**

Refraction Kit with Cylindrical Light Box

**Technical Specification :**

It is an apparatus for the path of light through 2 types of medium, namely water and air, with different densities. The light source is in water. Use 12 V AC/DC electricity.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003