

Product Name :
Refrigerator

Product Code :
NLAB-SCIENTIFICEQ12010



Description :

Refrigerator

Technical Specification :

Type: Chest

Temperature Range: between -18°C to -20°C

Capacity: 576 Liter

Class: HF 576 N

Power Connection: 240 V, 50/60 Hz

Provide the warranty period of 24 months

Storage of samples

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com