

Product Name :
Research Coaxial Microscopes

Product Code :
NLAB-EDUSCHOOLEQ110002



Description :

Research Coaxial Microscopes

Technical Specification :

Microscopes are latest for Research and Analytical work in Laboratories, Hospitals, Advanced Research Institutes and Colleges.

It is equipped with co-axial focusing system and halogen (optional LED) illumination for sharp and bright images.

The robust arm has co-axial focusing knobs on both sides which are based on a stage gear reduction system.

Travelling on ball bearing guides thus providing a smooth and precise motion control.

It has large base with built in 6V-20W Halogen bulb operates on 220V (optional 110V) and light source with ON/OFF switch with intensity regulator.

Large Mechanical Stage (145x125mm) fitted with low drive coaxial controls for precise x & y axis movements.

A movable sub stage condenser N.A.1.25 with iris diaphragm and filter holder.

Quadruple Nosepiece with rubber grip with following optical combinations: Achromatic or Semi Plan or Din Size Semi Plan Objectives: 4X,10X, 40X SL & 100X SL (oil).

Eyepieces: WF 10X Based on observation heads (interchangeable), it is available in the following versions.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003