

**Product Name :**  
Research Microscopes

**Product Code :**  
NLAB-EDUSCHOOLEQ170001



**Description :**

Research Microscopes

**Technical Specification :**

Research Microscopes are used for education, training analytical and research work in Colleges, Hospitals, Laboratories and institutes.

It has a heavy sturdy body & base for stability and is fitted with Quadruple revolving nosepiece.

Its stage size is 140x120mm with low coaxial controls for easy slide manipulation on both the axis with graduation and verniers, range of movements is 80x50mm.

Condenser unit is having rack & pinion movement with Abbe Condenser N.A.1.25, Iris-diaphragm and filter holder.

Having built-in base illumination system fitted with a Halogen Lamp 6V-20W operated on 220 volts (Optional 110volts).

Built in transformer is controlled by variable solid state electronic knob.

There are three different inclined observation heads (interchangeable): Monocular, Binocular & Trinocular inclined at an angle of 45° and rotatable through 360°.

Having following optical combinations: Achromatic or Semi Plan or Din Size Semi Plan Objectives: 4X,10X,40X

SL & 100X SL (oil) Eyepieces: WF 10X

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003