

Product Name : Resistance Coils with Sockets	Product Code : NLAB-EDUSCHOOLEQ1730204
--	--

--

Description :

Resistance Coils with Sockets

Technical Specification :

Two brass terminals are provided and the value of each coil is clearly marked on the top of the case.

The Coil is adjusted to accuracy with in 0.1% Any range from 0.1 ohms, 120 ohms.

Comprising a non-inductivity wounded coil of double silk covered constantan wire and enclosed in a circular moulded bakelite case.

Naugralabequipments
Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003