

Product Name : Resistance Coils with Sockets	Product Code : NLAB-EDUSCHOOLEQ1730204
Description :	
Resistance Coils with Sockets	
Technical Specification :	
Two brass terminals are provided and the value of each coil is clearly marked on the top of the case. The Coil is adjusted to accuracy with in 0.1% Any range from 0.1 ohms, 120 ohms. Comprising a non-inductivity wounded coil of double silk covered constantan wire and enclosed in a circular moulded bakelite case.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003