

Product Name :
Respironics Bipap Pro

Product Code :
NLABMEDICAL-M3390013

Description :

Respironics Bipap Pro

Technical Specification :

Compact, Lightweight and easily portable to fit well on patients' nightstands and with their lifestyles

Quiet so the patient and bed partner are not disturbed by noise

Ramp lowers the pressures at the start of therapy

Encore smart card records the date, time, ramp usage, duration of each use, and FOSQ answers

Advanced blower design and pressure feedback provide stable pressure and allow the system to tolerate most mask leaks

Compensates for altitude changes to provide comfortable, accurate pressure despite changes in altitude

Automatic on/off feature eliminates the need for a remote control yet still provides easy control over therapy

Mask removal indicator can be set to sound in case the mask is removed (may not function with all masks and circuit configurations or at high altitudes)

Dual time meters provide immediate status of total machine usage and cumulative patient compliance (time at pressure)

Universal power supply simplifies travel worldwide

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

