

Product Name :
Retention Tester

Product Code :
RET-IN-0001



Description :

Retention Tester

Technical Specification :

The is specifically designed for investigative testing of additive influences on pulp quality. The (retention tester) is an ideal and unique tool for controlling and setting up new chemicals, testing new wires, identifying additive problems and developing new paper products. This automated system offers a tight control on suspension shearing consistency, additive contact time optimization and drainage conditions. Because of this, the produces sheets that have technical characteristics closely related to those of industrial paper. A typical FRET test output yields a handsheet, the white water and the drainage vacuum curve for study.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773