

Product Name :
Retractable Fall Arrester Edge 15 Fall Arrester

Product Code :
NLAB-SAFETYLABP20056



Description :

Retractable Fall Arrester Edge 15 Fall Arrester

Technical Specification :

Wire rope type retractable fall arrester.

Maximum load capacity: 100 kg.

Strength: 25 kN.

15 metre in length.

Durable polyamide casing.

Included conduit rope.

4 mm galvanised steel rope / 4.8 mm stainless steel rope with swivel snap hook.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com