

**Product Name :**  
Retractable Fall Arrester Edge 20 Fall Arrester

**Product Code :**  
NLAB-SAFETYLABP20057



**Description :**

Retractable Fall Arrester Edge 20 Fall Arrester

**Technical Specification :**

Wire rope type retractable fall arrester.

4 mm galvanised steel rope / 4.8 mm stainless steel rope with swivel snap hook.

20 metre in length.

Heavy duty polyamide casing.

Strength: 25 kN.

Included conduit rope.

Maximum load capacity: 100 kg.

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com