

Product Name :
Retractable Fall Arrester Edge 2.5 Mtr Mini Block

Product Code :
NLAB-SAFETYLABP20053



Description :

Retractable Fall Arrester Edge 2.5 Mtr Mini Block

Technical Specification :

Fitted with an inertia breaking system with in-built shock absorber, snap hook.

2.5 metres in length.

Webbing type retractable fall arrester.

Weight:1 kg

47 mm wide webbing.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003