

**Product Name :**  
Retractable Fall Arrester Edge Compact

**Product Code :**  
NLAB-SAFETYLABP20054



**Description :**

Retractable Fall Arrester Edge Compact

**Technical Specification :**

30 mm wide webbing.

Swivel steel snap hook with fall indicator.

With shock absorber and snap hook.

Compact and lightweight.

Maximum load capacity 100 kg.

Webbing type retractable fall arrester.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

