

Product Name :
Retracting Microtome

Product Code :
NLAB-LABORATORYEM28023



Description :

Retracting Microtome

Technical Specification :

Accurate specimen advance.

Convenient specimen holder.

Protective cover for easy cleaning & lubrication.

Safety device.

Reduced hand stress.

Section thickness range: 1-50 microns,1-50 microns,1-20 microns

Section thickness setting: In steps 1 micron

Knife length: 120mm

Operation: Manual

Block (Specimen) Holder: Set of three

Object feed: With wheel

Vertical specimen stroke: 55mm,45mm

Max. specimen size: 32mm ,40x55mm

Adjustment with locks: Provided

Thickness control: Manual

Specimen orientation: 360°

Knife angle: 0 to 10°

Razor holder grip: Single,Double,Double

Dust cover: Provided

Instruction manual: Provided

L Moulds: Optional 50x25x20mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003