

Product Name : Riester Schiotz C Tonometer	Product Code : NLABMEDICAL-M3460004
Description :	
Riester Schiotz C Tonometer	
Technical Specification :	
High quality agate bearing for an extremely long service life.	
Precision measurement on a scale of 0 to 20 subdivisions and 0 to -1 sub-divisions.	
Perfect reading of the scale with red pointer.	
All vital parts are made of stainless steel, other parts are chrome-plated.	
Complete in a black deluxe case with velvet-look inserts.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003