

**Product Name :**  
Riester Schiotz C Tonometer

**Product Code :**  
NLABMEDICAL-M3460004

**Description :**

Riester Schiotz C Tonometer

**Technical Specification :**

High quality agate bearing for an extremely long service life.

Precision measurement on a scale of 0 to 20 subdivisions and 0 to -1 sub-divisions.

Perfect reading of the scale with red pointer.

All vital parts are made of stainless steel, other parts are chrome-plated.

Complete in a black deluxe case with velvet-look inserts.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003