

**Product Name :**  
Ripple Tank Kit

**Product Code :**  
NLAB-EDUSCHOOLEQ1730165

**Description :**

Ripple Tank Kit

**Technical Specification :**

For the visual demonstration of the properties of wave motion and the effects of refraction, diffraction and interference.

High quality ripple tank kit physics laboratory instrument.

Specification:

halogen lamp working voltage: 12V (100W)

working voltage: DC 12V  $\pm$ 5%

power supply: 220V/50Hz

Weight: 10.5kg

Size: 35 x 30 x 45 cm

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003