

Product Name :
Rocking Microtome

Product Code :
NLAB-LABORATORYEM28024



Description :

Rocking Microtome

Technical Specification :

Cambridge rocking microtome is very well known and one of the best established microtomes in the world, which has given outstanding service to numerous microtomists through the years. It is designed for sectioning comparatively soft, wax-embedded tissues with a thickness range of 2 mm to 24 mm, in steps of 2mm. Calibrated in 2 micron divisions, the section cutting thickness can be set easily and quickly by a single control. It is recommended for elementary biology & pathology students. The rocking microtome is an extremely reliable instrument comprising of three moving parts. The lead screw mechanism incorporates quality precision engineering and the cutting stroke is controlled by a single control lever which minimizes the operator skill required. It carries wide choice of specimen holders cylindrical holders, standard orientating & wood block holders for sectioning wax embedded specimens.

Overall Dimensions 390x200x200 mm Nett Weight – 4 kgs.

Size of the cabinet 450x240x260 mm Gross Weight – 10 kgs.

The instrument is supplied with the following in a storing box :

Three Standard Specimen Holders

Razor with built in handle in a box

Operating Instructions and Vinyl Dust Cover.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003