

Product Name :
Rockness Hardness Tester

Product Code :
NLAB-MECHANICALAB400002



Description :

Rockness Hardness Tester

Technical Specification :

Rockwell Hardness Testing Machine is designed for measuring hardness of metals & alloys of all kinds, hard or soft, whether round, flat or irregular in shapes.

This machine is ideally suitable for laboratories, tool rooms, Heat treatment shops, R&D departments, inspection departments, foundries & educational institutions.

Automatic weight selection with automatic zero setting dial gauge.

Rockwell test minor load is 10 kgf & major loads are 60, 100, 150 kgf.

Rockwell hardness scales such as HRA, HRB, HRC, etc. is obtained by using different types of indentors(Diamond / Ball).

Test height x Throat is – 215 x 132 mm

Extra test height & throat of 295 x 148 mm, Model : MRS-150 is also available. Motorized versions are also available.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003