

**Product Name :**  
ROTAP Sieve Shaker

**Product Code :**  
NLAB-TESTINGLABM109073



**Description :**

ROTAP Sieve Shaker

**Technical Specification :**

**ROTAP' Sieve Shaker**

To make process of sieving simpler and quicker, Rotap sieve shaker is useful , this produces circular shaking of the sieves. At the same time the sieves are tapped. The mechanism for imparting circular action and topping is oil immersed in a tank and in motorized. The shaker can accommodate upto 7 nos. Sieves of dia 200 or 300mm, 450mm Suitable operation of 230 V A.C. Single Phase.

**ACCESSORIES:**

Time switch adjustable from 0-60 minutes in 5 minute intervals.

**'ROTAP' Digital Sieve Shaker**

Same as above with Digital timer 0-99 x 1 minute.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003