

Product Name :
Rotary Microtome (Erma Type)

Product Code :
NLAB-LABORATORYEM28025



Description :

Rotary Microtome Erma Type

Technical Specification :

Heavy cast iron base with entire mechanism covered with a heavy cast iron protective dust cover. It incorporates a Ball & Flange type object holder for exact centering of specimen. Lateral Movement of Knife holder ensures use of complete knife edge. Microtome has a feed range of 1-25 microns in steps of 1 micron each. Supplied with Razor 120 mm long with back and Handle, Dust cover, Honing Plate, Object Holder (set of 3) and oil can. Razor supplied is made out of alloy steel & duly hardened for long lasting knife edge.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003