

Product Name :
Routine Teaching Microscope

Product Code :
NLAB-EDUSCHOOLEQ173044

Description :

Routine Teaching Microscope

Technical Specification :

Specifications:

Coarse and fine focusing

Monocular head inclined 45°, with 360° rotation

Built in 1.25 focusable Abbe condenser

x4, x10 & x40 DIN Achromatic objectives

Quadruple nosepiece

LED substage illumination

Iris diaphragm with filter holder and blue filter

x10 wide field eyepiece

Focusing Control: Coarse & Fine

Maximum Magnification: x400 (x1000)

Illumination: LED, variable

Condenser: Abbe

Diaphragm: Iris

Objectives: x4, x10, x40

Ocular Type: Monocular

Eyepieces: x10 W

Nosepiece: Quadruple

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003