

Product Name :
RTD Packed Bed Reactor Apparatus India

Product Code :
NLAB-VOCATIONALTVET29005



Description :

RTD Packed Bed Reactor Apparatus

Technical Specification :

Feed Tank: Material Stainless Steel, Capacity - 20 Ltrs.

Feed Circulation: By compressed air.

Flow Measurement: Rotameter.

Reactor Column: Material Borosilicate Glass

Packing: Rasching Rings, Material Borosilicate Glass, Size 8-10mm (approx)

The whole set-up is ingeniously designed and schematically arranged on a powder-coated rigid structure.

Control Panel:-

Standard make On/Off switch, Mains Indicator etc.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003