

**Product Name :**  
Rubber Block Apparatus

**Product Code :**  
NLAB-ENGINEERINGLB26004



**Description :**

Rubber Block Apparatus

**Technical Specification :**

Rubber Block Apparatus Lab Equipments Manufacturer, Suppliers & Exporters Description of Rubber Block Apparatus:

Rubber Resistance Apparatus Consists of a panel on which a block of Rubber as well as a console equipped with ab Dial are mounted. The comparator shows the degree of deformation of the block. A total shear load of 12 kg is permitted.

Experimental Ability of Rubber Block Apparatus:

Observation of the driving force and the Shear deformation on a rubber block.

Demonstration of the relationship between the Shear and shear deformation.

Determination of the stiffness modulus for the rubber

Determination of the variation in the deviation of the applied load

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003