

Product Name :
Rubber Stopper

Product Code :
NLAB-SCIENTIFICEQ31007



Description :

Rubber Stopper

Technical Specification :

Solid rubber stoppers; soft, superior quality rubber resistant to acid, alkali and ammonia.

Number of cork indicated on the face.

Sizes available are:-

size no	top mm	bottom mm	length mm
0	13	9	15
1	14	11	18
2	16	13	20
3	19	14	21
4	21	15	22
5	22	17	23

6	24	20	24
7	26	20	24
8	26	21	27
9	28	21	31
10	30	26	33
11	33	30	33
12	34	30	34
13	38	31	37
14	41	32	42
15	43	34	40
16	45	35	36
17	49	40	40
18	53	48	33
19	58	44	45
20	60	51	56

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003