

Product Name :
S. S. G. N. Pan

Product Code :
NLAB-VOCATIONALTVET270118

Description :

S. S. G. N. Pan

Technical Specification :

Size : 1/1x150

Capacity : 20 Litres

Color : silver

Material : S.S.

Usage : food pan

G. N. pan use all hotels, bakery and kitchens.

Mixing the waffle batter in G. N. Pan

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003