

Product Name :
Salt Bridge Tube

Product Code :
NLAB-SCIENTIFICEQ18008



Description :

Salt Bridge Tube

Technical Specification :

For use in experiments involving the measurement of cell EMF under conditions of forward, zero, and reverse current flow. The bridge comprises an inverted 'U' tube with its ends closed by sintered discs. This apparatus provides a more effective method of making reversible cells than the wet filter paper or agar gel methods. Arm diameter 13mm, width 8cm and overall height 10cm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003