

Product Name :
Sand Density Cone Apparatus

Product Code :
NLAB-TESTINGLABM109026



Description :

Sand Density Cone Apparatus

Technical Specification :

The Sand Density Cone and accessories are used for determination of the insitu (in-place) density of soils. The apparatus is restricted to tests in soils containing particles not larger than 2" (50.8 mm) in diameter. The Sand Density Cone consists of a double cone assembly with a brass valve and comes supplied with a one-gallon (4,000 cm³) capacity sand jug. The Field Density Plates are used in conjunction with the Sand Density Cones. The center opening of the plate can be used as a template to control the diameter of the hole and the plate itself is used as a shelf to prevent the loss of soil during the removal operation.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003