

Product Name :
Sand Equivalent Device

Product Code :
NLAB-EDUCATIONTVET13022



Description :

Sand Equivalent Device

Technical Specification :

Standards EN 933-8, NF XP18-598

Including:

Plexiglass cylinders calibrated at 100 and 380 mm

500 mm graduated steel ruler

Rubber stopper

Plastic funnel

200 ml calibrated box

5 l plastic bottle

Washer tube with plug and siphon assembly

Calibrated piston for sand level measurement

2 mm sieve, diam. 200 mm

Concentrated solution 1000 ml

Stopwatch

Stand 1000 mm

Transport box

Metal funnel

Electric stirrer for sand equivalent test, speed 175 strokes / min, amplitude 8 "(203 mm), with timer and automatic stop, 220 V / 50 Hz CE compliant power supply

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003