


**Product Name :**  
Sand Pouring Cylinder

**Product Code :**  
NLAB-TESTINGLABM109025



**Description :**

Sand Pouring Cylinder

**Technical Specification :**

The soil testing equipment is used for determination of the dry density of compact, fine, medium grained soils and for layers not exceeding 50cm thickness. A circular hole is dug into the ground, all the soil from within it collected, weighed and dried, and the hole back filled with a standard uniform sand or fine gravel, poured from a calibrated container for calculating the volume of hole. The complete apparatus consists of a Sand Pouring Cylinder fitted with Conical funnel and Shutter, Cylindrical Calibration Container, and a Metal tray with a central hole.

Height of Cylinder	Diameter of Container	Height of Container
-----------------------	--------------------------	------------------------

300mm	160x40mm	
-------	----------	--

300mm	160x45mm	
-------	----------	--

205450x50mm

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003