

Product Name :
SAYBOLT VISCOMETER APPARATUS

Product Code :
NLAB-ENGINEERINGLB16012



Description :

SAYBOLT VISCOMETER APPARATUS

Technical Specification :

DESCRIPTION:-

Used for determining the Viscosity of fuel oil and lubricating oil both thin & thick. The apparatus consists of Stainless Steel bath with oil cup which is centrally placed in a water bath. The bath has a lid which contains a water cooling tube, two handle with two stirrer blades, thermometer socket and a straight heater, stirring is done by turntable arrangement. A thermometer can also be inserted into the cup cover. The jets one universal and one Furol can be screwed to the cup by a handle (without strainer, withdrawal tube and glassware) can operate on 220 Volts AC mains.

SERVICE REQUIRED:-

- Fuel Oil OR Lubricating Oil sample.

- Suitable Water.
- Suitable Table space to mount apparatus.
- 15 Amps. Single Phase Electrical Connection

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003