

Product Name :
Scada System Heat Transfer Lab

Product Code :
NLAB-ELECTRICALAB170001



Description :

Scada System Heat Transfer Lab

Technical Specification :

Technical Description

Scada System Heat Transfer Lab

The Heat Transfer Model, designed as an accessory for the , is used for the experimental investigation of heat transfer from a heated rod to an air flow. It comprises a rod oven and a copper specimen. The copper rod is heated in the oven to approx. 110°C, before it is placed in the Computer Linked Air Flow Bench. A temperature sensor is built into the specimen; this is connected to the oven and from there transmits its measured data to the Data Logging Unit. By measuring the cooling curve in a flow of air, the coefficient of heat transfer can be determined.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003