

Product Name :
School Lab Alloy Boss Head

Product Code :
NLAB-EDUSCHOOLEQ11004



Description :

School Lab Alloy Boss Head

Technical Specification :

Make sure equipment is secured in place on a retort stand with this simple, effective clamp.

Nickel plated brass tommy screws.

Die-cast in the natural alloy.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003