

Product Name :
School Lab Double Burette Clamp

Product Code :
NLAB-EDUSCHOOLEQ11013



Description :

School Lab Double Burette Clamp

Technical Specification :

Use two burettes simultaneously with the double burette clamp.

Spring is positive, allowing easy operation without damaging glassware.

This strong clamp has plastic jaws and is sensitive enough not to damage glassware.

Suitable for burettes up to 100ml capacity.

High strength casting.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003