

Product Name :
School Lab Mini Bipod Stand

Product Code :
NLAB-EDUSCHOOLEQ11028



Description :

School Lab Mini Bipod Stand

Technical Specification :

The mini bipod stand is a stable platform that can be used to hold a beaker, funnel or watch glass.

57 x 57mm top.

70mm high.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003