

Product Name :
School Lab Polypropylene Beaker

Product Code :
NLAB-EDUSCHOOLEQ11040



Description :

School Lab Polypropylene Beaker

Technical Specification :

Avoid broken glass in the lab by using polypropylene beakers.

Polypropylene.

Graduated.

This beaker has a spout and graduation marks.

Spouted.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003