

Product Name :
School Lab PVC Tubing

Product Code :
NLAB-EDUSCHOOLEQ11044



Description :

School Lab PVC Tubing

Technical Specification :

The flexible, clear PVC tubing can be sterilized up to 125 degrees Celcius.

The 3mm works perfectly for pneumatics and hydraulics experiments with syringes.

Elastic and flexible.

Can be sterilized up to 125°C.

Non-toxic.

Length: 1m.

Available in the Following Internal Wall Thickness:

3.0mm 0.75mm

5.0mm 1.50mm

6.5mm 1.50mm
8.0mm 2.00mm
10.0mm 2.00mm
12.5mm 2.25mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003