

Product Name :
School Lab Solid Rubber Stopper

Product Code :
NLAB-EDUSCHOOLEQ11053



Description :

School Lab Solid Rubber Stopper

Technical Specification :

It is good to keep a variety of sizes and types of stoppers available for the various experiments taking

Size guide for solid rubber stopper:-

12.7 x 18.2 x 26.5mm (bottom x top x length)

21.6 x 25.5 x 24.5mm (bottom x top x length)

14 x 19 x 21mm (bottom x top x length)

21 x 26.5 x 27.5mm (bottom x top x length)

28 x 32.9 x 29.6mm (bottom x top x length)

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com